DS202 1Ph Dual Stage Rotary Vane Pump

Model G3850-80240

Manuale di istruzioni Bedienungshandbuch Notice de mode d'emploi User Manual

87-900-143-01 12/2013

Notices

© Agilent Technologies, Inc. 2011

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

Publication Number: 87-900-143-01

Edition

Edition 12/2013
Printed in ITALY

Agilent Technologies Italia S.p.A. Vacuum Products Division Via F.Ili Varian, 54 10040 Leinì (TO) ITALY

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or

contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Trademarks

Windows and MS Windows are U.S. registered trademarks of Microsoft Corporation.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

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
Original Instructions

General Information

This equipment is destined for use by professionals. The user should read this instruction manual and any other additional information supplied by Agilent before operating the equipment. Agilent will not be held responsible for any events occurring due to non-compliance, even partial, with these instructions, improper use by untrained persons, non-authorized interference with the equipment or any action contrary to that provided for by specific national standards.

The DS202 1Ph pump is a dual-stage, rotary vane pump oil sealed, driven by a single-phase electric motor.

This high vacuum pump is suitable for pumping non corrosive gases.

The following paragraphs contain all the information necessary to guarantee the safety of the operator when using the equipment. Detailed information is supplied in the appendix "Technical Information".

This manual uses the following standard protocol:

WARNING!

Warning messages are for attracting the attention of the operator to a particular procedure or practice which, if not followed correctly, could lead to serious injury.

CAUTION!

The caution messages are displayed before procedures which, if not followed, could cause damage to the equipment

NOTE

Notes contain important information taken from the text.

Symbols used

The following symbols are used consistently throughout in all illustration:

Suction flange

Filling up operating fluid

Exhaust flange

Warning "see installation/operating mode instruction"

Gas Ballast valve

Hot surface "Danger of burns if hot parts are touched"

4 Instructions for Use Storage

Storage

When transporting and storing the pump, the following environmental requirements should not be exceeded:

- temperature: from -20 $^{\circ}$ to +70 $^{\circ}$ C
- relative humidity: 0 95 % (non-condensing)

Preparation for Installation

The pump is supplied in a special protective packing. If this shows signs of damage which may have occurred during transport, contact your local sales office.

Total weight of the pack, including the pump, is approx. 32 Kg.

When unpacking the pump, be sure not to drop it and avoid any kind of sudden impact or shock vibration to it.

Do not dispose of the packing materials in an unauthorized manner. The material is $100\,\%$ recyclable and complies with EEC Directive 85/399.

NOTE

Normal exposure to the environment cannot damage the pump. Nevertheless, it is advisable to keep it closed until it is installed in the system, thus preventing any form of pollution by dust.


Figure 1

4 Instructions for Use Installation

Installation

Do not install or use the pump in an environment exposed to atmospheric agents (rain, snow, ice), dust, aggressive gases, or in explosive environments or those with a high fire risk. During operation, the following environmental conditions must be respected:

- temperature: from +12 °C to +40 °C
- relative humidity: 0 95 % (non-condensing)

WARNING!

Protect against short circuits and overload by installing on Agilent Device electrical main line an automatic circuit breaker of proper capacity (see table here below):

Tab. 1

P/N	200/240 Vac 50Hz	200/240 Vac 60Hz
G3850-80240	8A	8A

CAUTION!

Before starting the pump, fill up with lubricating oil as the pump is delivered empty.

WARNING!

Take out the protective caps on the suction and exhaust flanges before doing anything else. In the event of an accidental start-up, the air inside the pump could violently expel the protective caps and harm the operator.

WARNING!

During installation, pay maximum attention that the suction flange is connected to the vacuum chamber and the exhaust flange is not closed (see the following figure). The pump must not be used as a compressor.


Figure 2

WARNING!

A safety device (indicated by the red arrow in the figure below) is installed on the pump to guarantee that the pressure inside the oil vessel never exceeds the maximum permitted value of 1.7 bar (abs). This is a safety device and for this reason the cap must not be replaced, damaged or tampered with. Temporary removal (by unscrewing) is allowed to fill oil into the pump when needed. Users are responsible for screwing the device back in its original position after filling and before starting the pump.

Failure to respect these precautions may cause damage to the machine and injury to the operator.

A mixture of intake gas and oil will issue from the safety cap if the pressure in the oil vessel exceeds the allowed limit values (e.g. obstructed exhaust). The flow will continue until the pump is stopped and/or the cause of overpressure is removed (e.g. obstruction is removed from exhaust).

Figure 3

Check that your electrical mains voltage corresponds to that indicated on the motor's plate.

Connect the pump to the power supply.

4 Instructions for Use Use

Use

Before being put into service, in order to reach maximum vacuum, the pump must be left running for about an hour with the gas ballast valve open. This will eliminate any humidity from the oil.

There are no special procedures for switching the pump on; it needs only to be connected to the electric power by means of the bipolar switch (version with single-phase motor).

WARNING!

The pump is designed for operation with neutral or non-corrosive fluids. It is absolutely forbidden to use potentially explosive or flammable substances.

There are no special procedures for switching the pump off; it needs only to be disconnected from the electric power by means of the bipolar switch (version with single-phase motor).

Maintenance

Personnel responsible for pump operation and maintenance must be well-trained and must be aware of the accident prevention rules.

WARNING!

Death may result from contact with high voltages. Always take extreme care and observe the accident prevention regulations in force.

When machine is powered take care on account of moving parts and high voltages.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of $60\,^{\circ}\text{C}$.

4 Instructions for Use Maintenance

WARNING!

Always disconnect the power supply to the pump before starting maintenance work. Place a special warning signs over the power supply breaker switch: MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON. When finished, remove the safety warning.

WARNING!

Do not change the oil immediately after stopping the machine as the oil temperature may still be high.

CAUTION!

When replacing a part, take great care.

NOTE

Before returning the pump to the constructor for repairs the "Request for Return" sheet attached to this instruction manual must be filled-in and sent to the local sales office. A copy of the sheet must be inserted in the pump package before shipping.

If a pump is to be scrapped, it must be disposed of in accordance with the specific national standards.

Disposal

Meaning of the "WEEE" logo found in labels

The following symbol is applied in accordance with the EC WEEE (Waste Electrical and Electronic Equipment) Directive.

This symbol (valid only in countries of the European Community) indicates that the product it applies to must NOT be disposed of together with ordinary domestic or industrial waste but must be sent to a differentiated waste collection system.

The end user is therefore invited to contact the supplier of the device, whether the Parent Company or a retailer, to initiate the collection and disposal process after checking the contractual terms and conditions of sale.

4 Instructions for Use Disposal

Technical Information

Section I 67
Lubrication 69
Vacuum Seals 70
Gas Ballast Valve 71

Anti-Suckback Device 72

Technical Data 78
Dimensions 80

Safety Precautions 82

Transport and Installation 83

Preliminary Operations 85

Section II 86
Electrical motors 86

Single Phase Universal Motors 86

Single-Phase Universal Motors Voltage Setting 87

Connections to the Inlet and Exhaust Flanges 88

Starting and Running the Pump 90

Stopping the Pump 92

Original Instructions

5 **Technical Information Disposal**

Safety Rules 93 Warning Notes 94 96 Caution notes Maintenance Actions 97

99 Lubricants

Scheduled Maintenance Cards 104

Troubleshooting and Corrective Action 148

Spare Parts List 152

Section I

Technical Description

The DS202 1Ph pump is a dual-stage, rotary vane pump oil sealed, driven by a single-phase electric motor.


Figure 4

This vacuum pump is suitable for pumping non corrosive gases.

The main features are:

- all parts in direct contact with the fluid pumped are free of copper alloys;
- all materials are carefully selected to provide extended life;
- a high capacity gas ballast device allow to pump condensable vapors;
- all the parts composing the pump are fully replaceable due to the close machining tolerances and to the centering obtained by using reference pins;

5 Technical Information Technical Description

 due to its design features and low number of gaskets, the pump requires little maintenance, disassembly and reassemble are easy and require minimal time.

The pump works with force-feed lubrication, provided by an auxiliary gear-pump driven by the rotor of the vacuum pump itself. This ensures proper lubrication even when pressures are close to atmospheric.

The entire pump functional block is immersed in the oil contained in the casing. The oil guar3antees perfect sealing of the discharge valves, enters the pump to ensure lubrication and sealing of the parts inside, facilitates heat dissipation and reduces pump noise.

The pump is equipped with a special anti-suckback device which automatically isolates the vacuum system when the pump stops. This avoids rises in pressure or oil flow in the vacuum system while air is allowed back into the stator chambers.

The air entering the pump after the anti-suckback device has closed prevents the oil in the casing from filling the stator chambers.

The inner seals are achieved by the lubricating oil, thanks to the close machining tolerances.

There is only one oil seal on the rotor shaft, the one preventing oil from seeping out of the pump. On the DS202 1Ph pump, this seal acts on a surface-hardened area of the shaft guaranteeing long rotor life.

The vacuum pump is connected directly to the electric motor through a flexible coupling, so that motion is transmitted even in case of poor alignment.

A fan fitted on the joint and housed in the pump support produces a forced air flow over the finned surfaces of the casing to avoid oil overheating.

Lubrication

The pump's lubricating system is force-feed type.

Oil circulation is obtained by means of a gear pump connected to the rotor shaft.

The oil is drawn through a gauze filter to prevent any foreign bodies from entering the pump.

Shunt-mounted on the delivery line is a hydraulic piston which, besides actuating the anti-suckback device (see ANTI-SUCKBACK DEVICE), also regulates pressure by discharging excess oil flow directly into the tank.

The required flow of oil under pressure passes through ducts that are drilled in the walls of the pump and lubricates the bushings and the inside parts.

Vacuum Seals

A special feature of this pump is the low number of gaskets that are employed.

The seals in the circuit are obtained by means of VITON gaskets.

The careful surface finish of the various parts of the pump means that vacuum sealing of the functional block is ensured by the film of oil separating metal surfaces. Sealing of the rotor shaft is guaranteed by a rotating gasket with dust-guard lip.

The suction flange and duct are sealed by mean of OR gaskets.

Gas Ballast Valve

When the pump sucks in vapors, these condense during compression and mix with the oil, forming an emulsion (an aqueous vapor, for example) or a solution (organic solvent vapors, for example).

A number of problems arise from this, such as the impossibility of obtaining high vacuums, the alteration of the properties of the oil, could cause scaling on parts of the pump.

To avoid this, during compression at the second stage, atmospheric air is let into the pump through an adjustable valve, or "gas ballast valve", located at the top of the pump.

In this way, the discharge valve of the second stage opens through the effect of the atmospheric air before the partial pressure of the vapor reaches saturation point, thus preventing condensation from occurring. The vapors are expelled mixed with air.

Anti-Suckback Device

The pump is equipped with a special anti-suckback device to avoid air pressure rises and/or oil backflow towards the evacuated vessel when the pump is switched off. This device has a shutter which automatically closes the suction duct.

In this way the pump and vacuum system are completely isolated from each other and air can enter the pump without any risk for the vacuum produced in the system.

The device includes some special features, namely:

- drive obtained avoiding any form of contamination of the inlet duct by fluids (oil and/or air) used to command the shutter.
 Thanks to this, when the pump is started again, the pumpdown to vacuum conditions is extremely fast as these contaminants are not present and no degassing is therefore required;
- suction flange maintained closed even when the pump is idle, so that pollutants cannot enter from the environment and no oil can overflow from the pump.


Figure 5 View of the pump

5 Technical Information Anti-Suckback Device

Tab. 2

N°	DESCRIPTION	DS202	
1	Electric motor		
2	Rubber foot		
3	Screw		
4	Screw		
5	Screw		
7	Spacer		
8	Motor support		
9	Handle		
10	Screw		
11	Right sideplate		
12	Board		
14	Screw		
15	Support plate		
16	Cylindrical pin	(2) - (1)	
17	Screw		
18	Rubber crown	(2) - (1)	
19	Half-joint with fan		
20	Half-joint (pump side)		
21	Oil-seal ring	(2) - (1)	
22	Special screw		
23	OR gasket	(2) - (1)	
24	Gas Ballast Valve knob	SR03700219	
25	Half-plate		
26	Inlet screen	SR03700237	
27	OR gasket	(2) - (1)	
28	Suction flange		
29	Anti-suckback shutter	(2) - (1)	
30	Piston		
31	OR gasket	(2) - (1)	

N°	DESCRIPTION	D\$202	
32	Casing gasket	(2) - (1)	
33	Inner plate	SR03700260	
34	Pump support		
35	Anti-wear bushing		
36	OR gasket	(2) - (1)	
37	Kei motor side	SR8999980101	
38	First stage vane	(1)	
39	First stage rotor	SR03700319	
40	Spring	(2) - (1)	
41	Tongue	SR8999980001	
42	Second stage vane	(1)	
43	Second stage rotor	SR03700269	
44	Tie-rod	SR03700308	
45	OR gasket	(2) - (1)	
46	Bracket		
47	Screw		
49	Valve cover plate		
50	Valve spring	(2) - (1)	
51	Valve gasket	(2) - (1)	
52	First stage stator	SR03700316	
53	Middle plate	SR03700262	
56	GasBallast Valve shutter	(2) - (1)	
57	GasBallast Valve spring	(2) - (1)	
58	Second stage stator	SR03700321	
59	Oil pump joint	SR03700280	
60	Special screw		
61	End plate	SR03700285	
62	Tongue		
63	Washer		
64	Nut		
69A	External oil pump gear	SR03700281	

5 Technical Information Anti-Suckback Device

N°	DESCRIPTION	DS202	
69B	Internal oil pump gear	SR03700277	
70	Oil pump cover		
71	Screw		
73	Plate		
74	Filter		
75	Casing		
76	Exhaust flange		
77	OR gasket	(2) - (1)	
78	Oil fill plug	03.701218	
79	Oil drain plug	SR03700256	
80	Oil drain plug gasket	SR03700256	
81	Left sideplate		
83	Spring		
84	Filter		
85	Flange		
86	Sight glass	03.701326	
87	OR gasket	(2) - (1)	
99	Special screw		
100	Piston		
101	Arm		
102	Spring	(1)	
103	OR gasket	(2) - (1)	
104	OR gasket	(2) - (1)	
105	Sleeve		
106	Spring	(1)	
107	Ring	(2) - (1)	
108	Arm		
109	Insert		
110	Special washer		
111	Screw		

Tab. 3

1	Part of DS 202 Major Maintenance Kit (P/N9499381)
2	Part of DS 202 Minor Maintenance Kit (P/N 9499370)

Technical Data

The following table list the main technical data of the DS202 1Ph pump. $\,$

Tab. 4

TECHNICAL DATA	Hz	UNITS	DS202
FREE AIR DISPLACEMENT	60	I/min (cfm)	192 (6,8)
	50	I/min (m³/h)	160 (9,6)
PUMPING SPEED*	60	cfm	5,8
	50	m³/h	8,3
ULTIMATE PARTIAL PRESSURE*		mbar	10-4
ULTIMATE TOTAL PRESSURE*		mbar	2.10-3
ULTIMATE TOTAL PRESSURE WITH GAS BALLAST*		mbar	2.10-2
WATER VAPOR TOLERANCE		mbar	15
WATER VAPOR CAPACITY		g/h	100
OIL CAPACITY min/max		I	0,4/0,6
MOTOR RATING 1ph	50/60	kW	0,37/0,37
MOTOR RATING 3ph	50/60	kW	
NOMINAL ROTATION SPEED	50	rpm	1450
	60	rpm	1730
OIL TEMPERATURE		°C	50
(pump operating)**		°F	122
Installation category	II		
Pollution degree	2		
OPERATING TEMPERATURE RANGE		°C	12 - 40
WEIGHT		Kg	24
		lb	53
INLET FLANGE		DN	25KF
EXHAUST FLANGE		DN	25KF
Dimensions (cf. also next table):			

TECHNICAL DATA	Hz	UNITS	DS202
- length, 1ph		mm	491
- width		mm	134
- height		mm	230
Nominal voltages:	50	V	200-240
	60	V	200-240
Maximum currents:	50	А	8.0 2.7
	60	А	8.0 2.7

^{*} According to PNEUROP 6602

^{**} At ultimate total pressure, 20 °C (68 °F) room temperature

^{***} For more details about the electrical motors see pages from 17 to 20

Dimensions

Tab. 5

MODEL	UNITS	A	В	C	D	Н	L	M	N	P	R	w
DS202	mm inch											


Figure 6 Pump dimensions


Figure 7 Typical Performance Curves

Safety Precautions

- Always carry the pump by means of the handles provided.
- The pump must be set in position taking the upmost care in order to avoid accidental falls.

WARNING!

In case of a need to handle the pump after a period of operation, it must be left to cool first as the external surface temperature may be in excess of $60\,^{\circ}\text{C}$.

Transport and Installation

The pump is shipped to the customer inside cardboard boxes.

Total weight of the pack, including the pump, is 31 Kg.

The case must be handled with care, using appropriate lifting equipment.

CAUTION!

When moving the case, ensure that it is securely bound to the lifting equipment and that the equipment is strong enough to support the weight.

The pump's working environment is a traditional industrial environment. Naturally sites with corrosive vapors or excessive heat are best avoided.

Room temperature should ideally be between 12 °C and 40 °C.

If the temperature is not inside this range, consult Agilent technical service for the changes required.

5 Technical Information Transport and Installation

Setting the pump in position should be performed as follows:

- Pump laid on the ground. There are no special instructions for this type of installation, except that the floor should be as flat as possible and suited to bear the weight of the pump (it should ideally be a concrete floor) and of any accessories mounted on it. Note that the pump is stable on its base plate and it should not be necessary to anchor it to the floor with bolts and screws; also vibrations to and from the pump are greatly reduced by the use of rubber feet.
- Pump off the ground. In this case, the user must design a suitable support structure, remembering the following points:
 - the plane supporting the pump must be perfectly horizontal;
 - the structure should be adequately rigid;
 - the relevant safety precautions should be applied.
- Note also that the pump should be attached to the supporting structure after replacing the rubber feet with special antivibration feet, which should be screwed to the pump base and to the supporting plane.

After taking the pump out of its packing case, you are advised to make the following checks:

- **a** Ensure that the pump has not suffered any damage during shipping.
- **b** Check that the guards are mounted correctly (side plates, refs. 11 and 81 in fig. *View of the pump*) and that there are no uncovered or loose parts.

Preliminary Operations

Before starting the pump, fill up with lubricating oil as the pump is delivered empty.

NOTE

A tin of oil is included in the packing.

For details on how to fill up, see Scheduled Maintenance Card 01.

CAUTION!

Oil must be poured into the casing through the special threaded plughole and NOT through the suction line.

WARNING!

Take out the protective caps on the suction and exhaust flanges before doing anything else. In the event of an accidental start-up, the air inside the pump could violently expel the protective caps and harm the operator.

Section II

Electrical motors

Single Phase Universal Motors

They are world wide motors, dual voltages and dual frequencies, and they are in accordance with major international standards (UL, CSA, CE).

The motor has a main on/off bipolar switch.

The motor includes a thermal protection against overloads. If the motor is switched off by this relay, it can be started again after a few minutes.

Consult the trouble shooting section in order to understand and remove the overload cause.

NOTE

If fuses are installed on the power line, they must be in accordance with the table below.

If you start a pump at low temperature, the current absorbed by the motors will be for several seconds higher than the nominal one, for this reason the fuse on the line must be of the slow-blow type.

Tab. 6

PUMP	NOMINAL VOLTAGE [V]	HZ	RPM	kW	START UP Current [A]	FULL LOAD CURRENT [A]	MAXIMUM FUSE RATING [A]
Ds 202	200 to 240	50/60	1730	0.5	24	4.4	6

Single-Phase Universal Motors Voltage Setting

Before connecting to the mains, check that your electrical mains voltage corresponds to the motor voltage setting.

WARNING!

Disconnect the motor from the mains before opening the electrical box.

The motor can only operate at 200-240 V (50/60 Hz) nominal voltages.

WARNING!

Close the motor electrical box cover before connecting to the main.

Connections to the Inlet and Exhaust Flanges

Remove the protective caps from both flanges. Connect the system to be evacuated to the inlet flange, using a centering ring with OR and a locking collar.

NOTE

For guaranteed reliable sealing, use an OR gasket in Perbunan or Viton.

The inlet duct is equipped with a sieve filter preventing solid particles from entering and damaging the pump.

NOTE

When the gases to be pumped out contain dust, it is advisable to insert a dust filter before the inlet flange.

NOTE

When the gases to be pumped out contain large quantities of vapor, it is advisable to include a condense separator before the inlet flange.

To make best use of the pump's capacity, use only short, straight piping, with a diameter not smaller than that of the inlet flange.

NOTE

If rigid piping is used, it is good practice to use a flexible joint in order to avoid undue forcing of the connection on the pump.

The exhaust duct must be connected to a pipe that will take away the pumped out gases.

NOTE

Application of an oil trap filter is necessary to avoid pollution of the surrounding atmosphere by the oil present in the exhaust duct during pump operation.

CAUTION!

Never block the pump exhaust line. This would cause overpressure in the casing with the risk of breaking the glass window of the level indicator and/or expelling the oil seal gasket.

Starting and Running the Pump

Before being put into service, in order to reach maximum vacuum, the pump must be left running for about an hour with the gas ballast valve open. This will eliminate any humidity from the oil.

WARNING!

The pump is designed for operation with neutral or non-corrosive fluids. It is absolutely forbidden to use potentially explosive or flammable substances.

NOTE

If the pump is started with cold oil, initially more than normal noise will be heard; this will last for a few minutes only until the oil reaches its working temperature.

There are no special instructions for normal operation of the pump, which is delivered to you after completion of a running-in cycle in the factory. The maintenance operations for a pump in working condition may be found on the Scheduled Maintenance Cards.

If the pump is sucking in environments with liquids or vapors, it is important to leave it to run with the gas ballast open so as to avoid condensation of the vapors and contamination of the lubricating oil as a result. When there is no danger of the oil being contaminated, the gas ballast valve should be closed to reduce the pump's power absorption.

The gas ballast valve will prove useful in the cold season when, due to the lower room temperatures, it takes longer for the pump to reach rated temperature when switched on. Compression of the air forced in by the gas ballast valve develops heat, which in turn increases temperature of the oil faster.

NOTE

If gases without vapors are sucked in, you are advised to open the gas ballast device from time to time in order to eliminate the traces of humidity they contain.

NOTE

For repetitive work cycles, with brief time intervals in between, it is better not to stop the pump.

Stopping the Pump

There are no special procedures for switching the pump off; it needs only to be disconnected from the electric power by means of the bipolar switch.

When the pump is stopped, the anti-suckback device makes it possible to maintain vacuum in the vessel connected on the inlet flange of the pump.

If the pump is expected to be stopped for a lengthy period, or in any case if it has pumped in large amounts of vapors, it is good practice to run it with the gas ballast open and the inlet line closed for a few minutes before switching off in order to limit the risk of corrosion or scaling due to pollution of the oil by condensed vapors.

Safety Rules

Personnel responsible for pump operation and maintenance must be well-trained and must be aware of the accident prevention rules.

The accident prevention precautions contained in this section must be respected at all times during operation and maintenance of the pump to avoid damage to operators and to the pump.

These precautions are provided in the form of WARNING and CAUTION notes.

WARNING!

Operating procedures, technical information and precautions which, if not respected and/or implemented correctly may cause body harm to operators.

CAUTION!

Operating procedures, technical information and precautions which, if not respected and/or implemented correctly may cause damage to the pump.

Warning Notes

- **a** Death may result from contact with high voltages. Always take extreme care and observe the accident prevention regulations in force.
- b Always disconnect the power supply to the pump before maintenance work. Place a special warning signs over the power supply breaker switch: MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON.
- c If you are performing maintenance after the pump has been operating for a considerable time, allow sufficient time for it to cool as the external surface temperature may be in excess of 60 °C.
- **d** Failure to provide the pump with an earth connection may cause serious damage to operators. Always ensure that there is an earth connection and that it complies with the standards.
- When cleaning the pump and its component parts, avoid the use of flammable or toxic solvents, such as benzin, benzol, ether or alcohol. The recommendation is to use a soap and water solution, preferably in ultrasound washing machines, taking care to dry all the cleaned parts at temperatures under 100 °C in order to eliminate residual moisture.
- f Prolonged overloads or breakdowns may cause the electric motor to overheat, and to release noxious smoke; remove the power immediately as a precaution and do not approach the pump at least until you have provided ventilation to drive out the smoke. Take care not to breathe in the fumes remaining inside the pump in the course of repair work.
- **g** In case of fire, do not throw water on the pump. Switch the power off and use CO2 extinguishers.
- h Carefully inspect the flanges to ensure that there are no dust, oil, dirt or defects of the mating surfaces, before making the required connections.
- i Ensure that all joints and couplings are locked correctly before starting the pump again after repair work.

- j Do not wear any objects that may become entangled in the mechanisms and/or act as conductors (chains, bracelets, etc.).
- k Ensure that the tools to be used are in perfect working condition and have insulating grips, where necessary. Check that the insulating material of the cables and that the conductors of the test equipment do not show any signs of damage.
- I Do not replace the oil immediately after stopping the machine as the oil may still be at high temperature.
- m Perform repairs in clean and, where possible, dust-free areas. Protect all the clearances of connection points with suitable plastic caps and cover the machined surface areas of all parts stripped down until they are put back on the pump again.
- **n** Do not start the pump if the protective sideplates (refs. 11 and 81 of fig. *View of the pump*) are not properly in place.

Caution notes

- a Before putting the pump back into operation after a breakdown, inspect it and check carefully for any other signs of damage.
- b Use only tools that are in perfect working order and specially designed for the job; use of inappropriate or ineffective tools may cause serious damage.
- c Perform repairs in clean and, where possible, dust-free areas. Protect all the clearances of connection points with suitable plastic caps and cover the machined surface areas of all parts stripped down until they are put back on the pump again.
- **d** Always check the lubricant and that it is properly distributed through the pump; inadequate lubrication may damage the pump seriously.
- **e** Give the parts some form of marking as you strip them down to ensure that you reassemble them again in the proper order.
- f Check that there are no scratches or grooves on the machined shafts, in their seats inside the pump or on machine-ground surfaces. Slight scratches and abrasions may be eliminated with very fine emery paper or by a little light grinding.
- **g** Before putting a group together, always spread a little oil over inner parts and mating surfaces. Replace all seals with original spare parts before reassembling components.

Maintenance Actions

Maintenance may be seen as the totality of all scheduled and unscheduled maintenance work.

Maintenance includes:

- 1 SCHEDULED MAINTENANCE: Maintaining the nominal state of operation.
- 2 UNSCHEDULED MAINTENANCE: Restoring the nominal state of operation

NOTE

The frequency with which repairs are performed depends on the process and presence of substances that shorten pump life (dust, abrasives, solvents, water, chemically aggressive substances).

The pump must be cleaned at regular intervals of time.

CAUTION!

Do not clean with Alcohol the plastic or rubber components of the pump.

Use only the strictly necessary amount of lubricant; an excess of lubricating oil, like when there is none, may sometimes compromise proper operation of the pump.

Only the recommended lubricants, or lubricating oils with similar characteristics and known and experimented quality, should be used.

Oil changes must be made with the oil at a sufficiently high temperature, after leaving the pump to cool for a few minutes following operation.

The drain and filler plugs must not be left open any longer than is strictly necessary.

5 Technical Information Maintenance Actions

When performing maintenance, look out for all signals that may precede a breakdown, in particular:

- traces of corrosion;
- oil leaks;
- slack joints or couplings.

Maintenance technicians must:

- be aware of all applicable national directives concerning accident prevention during work on motor-driven pump and should know how to apply them;
- have read and understood all the sections on "Safety Rules";
- be familiar with the essential design features and operation of the pump;
- know how to use and consult the pump documentation;
- be concerned about proper operation of the pump;
- make a note of any irregularities in operation of the pump and take the necessary action, where appropriate.

Use original spare parts wherever possible and repair a broken part as best as possible on site or send it back to the manufacturer for repairs.

For all problems arising, or to order spare parts, refer to our service department.

Agilent Technologies Italia S.p.A. Vacuum Products Division Via F.lli Varian 54, 10040 Leini, (Torino) - Italy

Tel.: +39 011 997 9111 Fax: +39 011 997 9350

Toll-Free: 00 800 234 234 00

Lubricants

It will be readily understood how important adequate lubrication is to high technology pump like the Agilent vacuum pump. Correct use of appropriate lubricants makes a significant contribution to achieving best performance and warding off defects.

When handling lubricants, the following sanitary protection measures should be observed at all times:

- Avoid prolonged, excessive or repeated contact of the skin with products for lubrication, and also avoid directly inhaling the fumes or vapors of such products.
- Protect the skin by wearing appropriate clothes and equipment (e.g. special suits, glasses or, where permitted by the safety regulations, gloves) or by applying a special protective product.
- Clean the skin carefully after contact with the lubricants by washing freely with water and soap.
- Apply a skin cream after washing.
- Take off and change clothes or shoes on which oil has been spilled.
- Never put rags dripping with oil into the pockets of your clothes.

When disposing of waste lubricants, observe the following environment protection regulations:

- The lubricants risk contaminating the water and the ground! Therefore never pour lubricating products on to the ground, into water or in the sewage system. All violations of these rules are liable to persecution as provided for by law. When using lubricants always keep an oil can nearby.
- Take care in draining off waste oils. In disposal of these products respect all regulations in force concerning waste oil disposal.

The recommended lubricating oil is the Agilent Rotary Vane Fluid DS45 Type. The Rotary Vane Fluid DS45 Type is a general purpose mechanical pump fluid specifically engineered to provide superior performance in high speed direct drive mechanical pumps.

5 Technical Information Lubricants

These precisely distilled fluids (100 % solvent refined neutral paraffinic oil) deliver lower base pressure capability, faster pumpdown cycles, and reduced maintenance requirements on both the pump and the fluid.

It is absolutely necessary to continue using the lubricants initially used to fill the tank. If this is not possible for organizational or business reasons, use only products with the same characteristics as the previous oils.

Only use of lubricants of suitable quality will guarantee safe operation of the pump.

CAUTION!

Mineral oils and the PFPE oil are incompatible. To change from one type to another, the pump must be stripped down completely and all parts washed carefully to eliminate all oil residues.

If you expect to have to use other lubricants, first find out if the two products are compatible. In cases of doubt, the lubricant used up to that time must be flushed out by way of a pump flushing procedure.

CAUTION!

To avoid the risk of contaminating the oil, absolute cleanliness of the pump and surrounding area must be ensured during the lubrication procedures.

Tab. 7

PROPERTY	UNIT OF MEASURE	ROTARY VANE FLUID DS45 Type
Vapor pressure @ 25 °C	Torr	<1.10-7
Boiling point @ 0,01 torr	°C (°F)	133 (271,4)
Viscosity @ 40 °C	cSt (SUS)	53 (250)
Viscosity @ 100 °C	cSt (SUS)	7 (49)
Pour point	°C (°F)	-59 (-74,2)
Flash point	°C (°F)	264 (507,2)
Fire point	°C (°F)	288 (550,4)
Density @ 25 °C	g/ml	0,83

5 Technical Information Lubricants


Figure 8 Agilent Rotary Vane Fluid DS45 Type viscosity curve


Figure 9 Agilent Rotary Vane Fluid DS45 Type vapor pressure curve

Scheduled Maintenance Cards

Tab. 8

No.	FREQUENCY	PUMP STATUS (*)	OPERATION
01	MONTHLY	R	OIL LEVEL CHECKING
02	TWICE YEARLY	S	CLEANING THE INLET FILTER

(*) R=RUNNING S=STOPPED

NOTE

When the pump is equipped with:

- a filter for dust on the inlet side
- a condensation separator on the inlet side

these items must also be maintained, by following the Manufacturer's instructions.

NOTE

The numbers used in the drawings on the Scheduled Maintenance Cards are those used in Figure View of the pump.

		Scheduled Maintenance	Card No. 01
Pump type:	DS202 1Ph		
Frequency: Month	hly	Time: 5 min	Machine status: S R
Type of action:		LUBRICATION	
Item:		LUBRICATING OIL 949-9390M001	


Figure 10

No.	Description	Time	Equipment	Materials
1	Check the oil level and top up if	5 min	 Cleaning rags 	Oil P/N 949-9390M001
	necessary.		 Liquid detergent 	

WARNING!

Machine powered: take care on account of moving parts and high voltages.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.

 Check the oil level through the sight glass: there are two clearly visible arrows indicating maximum and minimum levels for proper pump operation.

NOTE

The amount of oil consumed depends on the type of work, particularly on the volume being evacuated and on the frequency of the vacuum cycles.

2. If necessary, open the pump filling plug on the top plane of the pump and add oil.

NOTE

Before starting the procedure, check that the quality and type of oil used corresponds to the indications of Table Agilent Rotary Vane Fluid DS45 Type characteristics of this manual.

CAUTION!

Take great care to avoid oily residues or dust from entering the pump and jeopardizing proper operation.

3. Close the fill plug and, using cotton or other suitable material rags, clean thoroughly the outer surface of the top plane of the pump.

CAUTION!

Never mix lubricants of different brands or quality. To change from one oil to another, the pump must be stripped down and all parts washed carefully, taking care to eliminate oil residues.

4. Open the gas ballast valve to purify the oil.

5 Technical Information Scheduled Maintenance Cards

		Scheduled Maintenance	Card No. 02	
Pump type:	DS202 1Ph			
Frequency: Twice Yearly		Time: 10 min	Machine status: S R	
Type of action:		CLEANING		
Item:		INLET FILTER		


Figure 11

No.	Description	Time	Equipment		Materials
1	Check the oil level and top up if	10 min	•	Cleaning rags	Oil P/N 949-9390M001
	necessary.		•	Liquid detergent	

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.

- 1. Remove the filter (26) from the inlet flange (28).
- 1. Clean the filter using an aqueous solution and dry with hot air.
- **2.** Put the filter back in place.

CAUTION!

The filter element is designed to prevent solid particles in the air from entering the inside of the pump and if it is not working properly, serious damage could be caused the pump's inner parts. In such cases, it may be worthwhile installing a large surface area filter in the pump's suction line to improve effectiveness and reduce maintenance frequency.

If the pump sucks in gases and vapors that contain large quantities of solid particles, the filter must be checked and cleaned more often than indicated here.

Tab. 9

UNSCHEDULED MAINTENANCE CARDS

No.	SUBJECT	FREQUENCY
101	General Warnings	
102	Changing the oil	Whenever necessary
103	Disassembly of oil level indicator	When damaged
104	Disassembly of anti-suckback device shutter	When damaged
105	Disassembly of electric motor and joint	When damaged
106	Disassembly of discharge valve	When damaged
107	Disassembly of oil circulating pump	When damaged
108	Disassembly of anti-suckback valve	When damaged
109	Disassembly of pump body	When damaged
110	Replacement of oil seal ring/wear bushing	When damaged
111	Disassembly of functional block rotors	When damaged

NOTE

The numbers used in the drawings on the Unscheduled Maintenance Cards are those used in Figure View of the pump.

		Unscheduled Maintenance	Card No. 101
Pump type:	DS202 1Ph		
Frequency: Twic	e Yearly	Time: 10 min	Machine status:
Type of action:		GENERAL WARNINGS	
Item:			


Figure 12

1. Pump disassembly and reassembly require some basic operations and use of tools found in most workshops.

NOTE

The pump's component parts are modular and therefore fully interchangeable. Whenever parts have to replaced, please specify in the order:

Pump type

Serial number (marked on pump identification plate).

Reference number on the detailed diagram, name and part number as indicated on the spare parts and accessories list (section SPARE PARTS LIST of this manual).

To simplify maintenance work, two kits have been devised.

The <u>Gaskets</u> Kit, containing all the gaskets recommended for substitution after disassembly of the pump and the most commonly used removable parts.

You are advised to keep this kit in a dry room, away from heat and light (sun and UV rays) in order to prevent hardening of the elastomer parts.

The <u>Maintenance</u> Kit, including, as well as the Gaskets Kit, the 1st and 2nd stage vanes and other parts that it may be necessary to replace when overhauling the pump completely. These parts have a life of approx. two years under conditions of normal operation.

Tab. 10

KIT	DS202 1Ph	ORDER CODE
GASKETS KIT	X	949-9370
MAINTENANCE KIT	X	949-9381

5. Before reassembling the pump, all parts must be cleaned thoroughly with a soap and water solution and rinsed with demineralized water.

Particular care must be taken in cleaning the lubrication ducts and oil passage nozzles, which must always be kept free of impurities.

CAUTION!

Take care that the soap and water solution does not remain in the various cavities and recesses of the pump, in the blind holes in particular.

- **6.** Then dry the parts, by heating where possible but without using compressed air (it generally contains polluting substances), taking care not to exceed a temperature of 80 °C.
- 7. The sequence of reassembly is the inverse of the disassembly operation. Take the same precautions:
- Spread some clean oil of the same type as used to lubricate the pump over the different parts, in particular the inner lip of the oil seal, the bushings and vanes.
- If the oil seal has the dust-block lip, 2/3 fill the area between the two lips with ISOFLEX TOPAS NB5051 grease.

CAUTION!

Do not allow the grease into the area of the circumferential spring.

- Lock without forcing the nuts on the four tie-rods (44) connecting the parts of the pump block. Before locking, check that all parts are touching as required.
- Before putting the springs and shutters of the discharge valves back in, pour a little oil into the stators through the holes in the valves.

CAUTION!

Before hooking the recall spring to the anti-suckback command bar, make sure you have put the bar all the way into the sleeve of the rotary feedthrough to avoid bending the tapered part of the bar.

NOTE

On DS202 1Ph pump, fit the joint on the motor shaft after careful removal of grease from the surfaces and application of LOCTITE 242.

		Unscheduled Maintenance	Card No. 102
Pump type:	DS202 1Ph		
Frequency: AS R	EQUIRED	Time: 30 min	Machine status: S R
Type of action:		CHANGING OIL	
Item:		LUBRICANT 949-9390M001	


Figure 13

Ma	Materials required:		Equipment required:	
•	Agilent Rotary Vane Fluid DS45	•	Fixed wrench for oil exhaust.	
	Type (1 liter) P/N 949-9390M001	•	Tank to collect oil.	
		•	Cotton (or other suitable material) rags for cleaning.	
		•	Soap and water solution for cleaning.	

Do not change the oil immediately after stopping the machine as the oil temperature may still be high.

WARNING!

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 14

NOTE

The oil has to be replaced:

- 1. when it has become particularly dirty due to the presence of contaminants in the environments being evacuated;
- 2. when the degree of vacuum achieved by the pump is no longer satisfactory, even after running for one hour with the gas ballast device open.
- 1. Place a tank for the waste oil under the oil drain plug (79).

NOTE

You are advised to empty the pump before it cools down completely so that the oil flows out more easily.

NOTE

Make sure the tank you use for the waste oil can hold up to 1 liter to avoid spill out.

- 8. Open the oil drain plug (79) slowly; let the oil flow out by opening the fill plug (78) located at the top of the pump. In the meantime, clean the inside of the drain plug (magnetic plug).
- **9.** After closing the drainage hole with its plug, make the pump perform a number of turns, giving it a brief current pulse, so as to let residual oil in the pump chambers flow out. Then drain off the residual oil again.
- **10.** After all the oil has been drained, close the drain plug tight; take the tank of waste oil away and clean thoroughly, using rags of cotton or other suitable material to dry.
- **11.** Clean thoroughly the pourer lip of the new oil tin; pour the oil slowly into the pump through the filler plug (78).

NOTE

Monitor through the sight glass at the front of the pump that you are not exceeding the maximum oil level permitted.

CAUTION!

For disposal of waste oil, it is strictly necessary to respect the currently applicable legislation on this matter.

12. To achieve a good vacuum level again, run the pump for one hour with the gas ballast open.

		Unscheduled Maintenance	Card No. 103
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 30 min	Machine status: S R
Type of action:		DISASSEMBLY	
Item:		SIGHT GLASS	


Figure 15

Materials required:	Equipment required:
OR gaskets (Gaskets kit p/n 949-9370)	 Mechanical tools
 Sight glass 	

5

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 16

2. Drain off the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

- **3.** Unscrew the securing screws (71) and disconnect the flange (85) from the casing.
- **4.** Take out the OR gasket (87) and the sight glass (86).

NOTE

When reassembling, replace the gaskets that you have stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil again following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 104
Pump type:	DS202 1Ph		
Frequency: WHEN DAMAGED		Time: 30 min	Machine status: S R
Type of action:		REPLACEMENT	
Item:		ANTI-SUCKBACK DEVICE SHUTTER	


Figure 17

Materials required:	Equipment required:	
■ Shutter	 Mechanical tools 	

WARNING!

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 18

- 1. Disconnect the flange (28) from the inlet pipe line.
- **13.** Remove the 6 securing screws (17) and remove the two halfplates (25) pressing the suction flange (28) into its seat.
- **14.** Remove the flange (28) and take out the piston (30) with the shutter (29) of the anti-suckback device, disengage it by exerting force on arm (101). Separate the shutter (29) from the piston rod (30) simply by pulling.

CAUTION!

Be careful not to scratch the sealing surface of the shutter (29).

15. Reassemble by proceeding in the inverse order.

		Unscheduled Maintenance	Card No. 105
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 30 min	Machine status: S R
Type of action:		REPLACEMENT	
Item:		ELECTRIC MOTOR AND JOINT	


Figure 19

Materials required:	Equipment required:
Electric motor	Mechanical tools

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 °C.


Figure 20

- 1. Disconnect the motor bearing plate (8) from the support (34) by removing the four screws (4).
- 2. Remove the half-joint/fan (19) from the support, together with the motor (1), the plate (8).

- **3.** 16. Slacken the nut blocking the half-joint (19) to the motor shaft and take it out, using leverage between the plate (8) and the back of the joint.
- **4.** Remove the four screws (5) and separate the plate (8) from the motor (1).
- **5.** Reassemble by proceeding in the inverse order.


Figure 21 Extractor

		Unscheduled Maintenance	Card No. 106
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 1 hour	Machine status: S R
Type of action:		DISASSEMBLY	
Item:		DISCHARGE VALVES	


Figure 22

Materials required:	Equipment required:
 Valves (Gaskets kit) 	 Mechanical tools

WARNING!

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 23

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

- 2. Disconnect the discharge flange from its piping.
- **3.** Remove the nuts (3) and remove the pump casing (75).
- **4.** Remove the securing screw (47) holding the splash guard (49) and remove it.
- **5.** Remove the two screws (47) on the valve press-bracket (46) and remove the bracket (46), the springs (50) and the shutters (51).
- **6.** Reassemble by proceeding in the inverse order.

NOTE

When reassembling, replace all the gaskets stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 107
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 1 hour	Machine status: S R
Type of action:		DISASSEMBLY	
Item: OIL CIRCULATION PUMP 949-9390M001			


Figure 24

Materials required:	Equipment required:	
Oil circulating pump P/N 949-9390M001	 Mechanical tools 	

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 °C.


Figure 25

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

- 2. Disconnect the discharge flange from its piping.
- **3.** Remove the nuts (3) and remove the casing from the pump (75).
- **4.** Remove the three screws (71) connecting the closing plate (70) to the side and take it out.
- **5.** Remove the pump gear-wheels (69).

NOTE

A magnet may be used to remove the gear-wheels.

6. Reassemble by proceeding in the inverse order.

NOTE

When reassembling, replace all the gaskets you have stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil again following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 108
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 1 hour	Machine status: S R
Type of action:		DISASSEMBLY	
Item:		ANTI-SUCKBACK DEVICE	


Figure 26

Materials required:	Equipment required:
Gasket kit	 Mechanical tools

WARNING!

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 27

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

- 2. Disconnect the discharge flanges from their piping.
- **3.** Remove the locking nuts (3) and remove the casing (75) from the pump.

- **4.** Remove the securing screw (47) and remove the splash guard (49).
- **5.** Remove the shutter from the anti-suckback device.

NOTE

For more details on how to remove the anti-suckback device, see Unscheduled Maintenance Card No. 104.

- **6.** Remove the screw (102), take out the arm (101) from the rod.
- 7. Release the spring (106) from the rod arm (108).
- **8.** Pull the rod (108) out of the rotary feedthrough (105) sleeve.
- **9.** Remove the two OR gaskets (104) from the rod (108) using a screwdriver, if necessary. Be careful not to bend the rod in doing so.
- **10.** Remove the sleeve (105) and remove it from the support (34) together with its OR gasket (103).
- 11. Remove the hydraulic piston (100) from its seat in the body.
- **12.** Reassemble by proceeding in the inverse order.

NOTE

When reassembling, replace all the gaskets stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil again following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 109
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 1 hour	Machine status: S R
Type of action:		DISASSEMBLY	
Item:		PUMP BODY	


Figure 28

Materials required:	Equipment required:
■ Gasket kit	 Mechanical tools

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 29

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

- **2.** Disconnect the discharge flanges from their piping.
- **3.** Disconnect the electric motor (1) complete with half-joint.

NOTE

For a description of how to disassemble the electric motor, see the Unscheduled Maintenance Card No. 105.

- **4.** Slacken the nut locking the half-joint (20) to the shaft of the rotor (39) and take out the half-joint. Remove the tongue (37) with it as well.
- **5.** Unscrew the securing screw (47) and remove the valve screen (49).
- 6. Slacken the four nuts (64) and take out the tie-rods (44) and the complete functional block (76), lifting it from the pump support (34) and take out the OR gaskets (27) and (45), from the gas ballast line and suction line respectively.

NOTE

Be careful not to cut the lip of the oil seal with the sharp edges of the shaft key seat.

7. Remove the screw (102), take out the arm (101) from the rod.

NOTE

When reassembling, replace all the gaskets stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil again following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 110
Pump type:	DS202 1Ph		
Frequency: WHEN DAMAGED		Time: 1 hour and 15 min.	Machine status: S R
Type of action:		REPLACEMENT	
Item:	OIL SHAFT SEAL / ANTI-WEAR BUSHING		BUSHING


Figure 30

Materials required:	Equipment required:
Oil shaft seal	 Mechanical tools
 Anti-wear bushing 	
 Gaskets kit 	

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

WARNING!

If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 $^{\circ}$ C.


Figure 31

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

2. Remove the worn out oil shaft seal (21) from the pump support (34) using a screwdriver for leverage and taking care not to scratch the seat.

NOTE

For a description of how to separate the functional block from the support, see the Unscheduled Maintenance Card No. 109.

- **3.** Remove the worn out oil shaft seal (21) from the pump support (34) using a screwdriver for leverage and taking care not to scratch the seat.
- **4.** Take the anti-wear bushing (35) off the shaft and take out the OR Gasket (36).

NOTE

For assembly of the new oil shaft seal, a special tool is needed.


Figure 32

- After applying silicone grease to the external diameter of the oil seal, engage it on the tool and push it into its seat.
- Reassemble by proceeding in the inverse order.

NOTE

When reassembling, replace all the gaskets stripped down.

CAUTION!

When reassembly is complete, proceed to fill the pump with oil again following the instructions of Unscheduled Maintenance Card No. 102.

		Unscheduled Maintenance	Card No. 111
Pump type:	DS202 1Ph		
Frequency: WHE	N DAMAGED	Time: 2 hours	Machine status: S R
Type of action:		DISASSEMBLY	
Item:		FUNCTIONAL BLOCK ROTORS	


Figure 33

Ma	terials required:	Eq	ipment required:
•	Oil shaft seal ring	•	Mechanical tools
•	Anti-wear bushing		
•	OR gaskets		
•	1st stage vanes		
•	2nd stage vanes		
•	Gas Ballast shutter		
	(the above items of material are included		
	in the Major Maintenance kit)		

WARNING!

Machine safety: disconnect all power lines. Hang the card "MACHINE UNDERGOING MAINTENANCE - DO NOT POWER ON" from the main breaker switch. When finished, remove the safety warning.

WARNING! If you have to perform maintenance on the pump after a considerable time in operation, leave it to cool as temperature of the outer surface may be in excess of 60 °C.

5 Technical Information Scheduled Maintenance Cards


Figure 34

1. Drain off all the oil from the pump.

NOTE

For a description of how to drain off the oil, see the Unscheduled Maintenance Card No. 102.

2. Detach the entire functional block from the pump support.

NOTE

For a description of how to separate the functional block from the support, see the Unscheduled Maintenance Card No. 109.

3. Disassemble the discharge valves.

NOTE

For a description of how to disassemble the discharge valves, see the Unscheduled Maintenance Card No. 106.

- **4.** Detach the end plate (61) from the rest of the pump body, exerting leverage if necessary on the reference notches between the parts. This will give you access to the retaining valve (56) and associated spring (57) at the gas ballast inlet.
- **5.** Disconnect the 2nd stage stator (58), remove the two vanes (42), take the 2nd stage rotor (43) off the shaft and take out the key (41).
- **6.** Remove the partition wall (53) from between the 1st and 2nd stages.
- 7. Detach the 1st stage stator (52) from the inner side (33), pressing back on the springs (40) of the vanes (38) to avoid the vanes from touching against the sharp corners of the body and becoming damaged.
- **8.** Take out the 1st stage rotor (39), the vanes (38) and the springs (40).
- **9.** Reassemble by proceeding in the inverse order.

NOTE

When reassembling, replace all the gaskets stripped down.

CAUTION!

After reassembly is complete, proceed to fill the pump again with oil following the instructions of Unscheduled Maintenance Card No. 102.

Troubleshooting and Corrective Action

NOTE

If the pump breaks down or does not work well, first check thoroughly that all the standards rules listed in the sections relevant installation and operation are currently observed.

If you cannot get sufficiently low pressure, the problem may lie with the system connected to the pump, i.e. the vacuum chamber and inlet line.

Even when the design sizes are correct, the problems listed below are still fairly frequent:

- a) faulty seals: check that the seals are fitted, are in good working orderand are clean; check all joints, taps, sight glasses, weld joints, etc.
- b) unsuitable materials: avoid the use of rubber, plastic or even metal parts that are porous or permeable.
- c) volatile substances: as far as possible, eliminate humidity, grease, dirt, rust, solvents (even small traces).

Tab. 11

	SYMPT0M	PROBABLE CAUSE	SUGGESTED REMEDY	CARD REF.
1	The pump does not start (or stops immediately).	A) The motor-pump coupling is defective.	Check state of the elastic joint and replace the parts as needed.	105
		B) Viscosity of the oil has increased.	Examine the oil and change if necessary (find out why it has deteriorated).	102
		C) Excessive friction inside pump.	Test the resistance to rotation, by turning the joint by hand. If the rotors are blocked, disassemble the pump body and check state of the parts inside.	111
		D) Electrical protections not set.	Adjust the protections.	
		E) Oil too viscous because of low room	Heat the oil up to a temperature of 20 °C or use	102

	SYMPTOM	PROBABLE CAUSE	SUGGESTED REMEDY	CARD REF.
		temperature.	an oil recommended for low temperatures.	
2	The pump stops (stalls).	A) Seizure of some internal parts.	Proceed as at 1-C.	
		B) Oil very dirty.	Proceed as at 1-B.	
3	The pump overheats.	A) Too much friction inside.	Proceed as at 1-C.	
		B) Oil ducts clogged inside.	Try changing the oil, otherwise disassemble the pump and clean it well.	102 111
4	Too much oil drained off by the exhaust.	A) Gas Ballast device broken (always open).	Check the device to ensure the gas ballast is properly adjusted.	111
		B) Heavy duty work (high suction pressures or very frequent cycles).	Check condition of the oil trap filter.	
		C) Leak in system.	Check the system (see "Note").	
5	The pump is loosing oil.	A) Seal on shaft defective.	Replace the oil seal and the anti-wear bushing.	110
		B) Other gaskets defective.	Identify the areas where seeping occurs and check the associated gaskets.	Refer to the card related to involved item
		C) Oil level window broken	Replace the oil level window.	

5 Technical Information Troubleshooting and Corrective Action

	SYMPT0M	PROBABLE CAUSE	SUGGESTED REMEDY	CARD REF.
6	The pump is very noisy even at steady operating conditions.	A) Exhaust valves damaged.	Inspect the valves.	106
		B) High suction pressure, or gas ballast device is open.	Check the system and the gas ballast valve (see items 4-A and 4-B).	
		C) Elastic star of joint is damaged.	Check and replace if necessary.	
		D) A screw of the joint has become loose.	Lock the screw after application of Loctite 222 (for slack threads).	
7	When the pump is stopped, the inlet pressure rises rapidly.	A) The anti-suckback device is defective.	Disassemble and check the anti-suckback device.	108
		B) The cause of the problem is not with the pump.	See "Note".	
8	The pump does not produce a sufficient vacuum.	A) Not enough oil, or oil contaminated.	Check level and appearance of the oil. Open the Gas Ballast for half an hour at least; change the oil whenever necessary.	01 102
		B) Oil ducts clogged inside.	Proceed as at 3-B.	
		C) Direction of rotation inverted.	Proceed as at 3-C.	
		D) OR gaskets communicating with the inlet line defective.	Replace the gaskets.	109
		E) Exhaust valves damaged or out of place.	Proceed as at 6-A.	

SYMPT0M	PROBABLE CAUSE	SUGGESTED REMEDY	CARD REF.
	F) Insufficient oil pressure in the lubricating circuit, preventing the antisuckback device from opening.	Inspect the gears of the oil pump and check that it is working.	107
	G) Internal moving parts worn.	Disassemble the pump, and replace any worn parts.	111
	H) Pump working with gas ballast device open.	Close the gas ballast device.	
	I) Anti-suckback valve fails to open.	Disassemble the suction flange and check the device.	104 108
	J) Parts of the functional block not in perfect contact because of deposits of dirt or failure to put a plug on the service holes.	Disassemble and check condition of the part surfaces and that plugs have been used.	111
	K) The cause of the problem is not with the pump.	See "Note".	

5 Technical Information Spare Parts List

Spare Parts List

Tab. 12

OTHER SPARE PARTS	0.TY for DS202 1Ph	P/N
MINOR MAINTENANCE KIT (*)	X	949-9370
MAJOR MAINTENANCE KIT (**)	X	949-9381

Tab. 13

ACCESSORY ITEM ON EXHAUST	0.TY for DS202 1Ph	P/N
OIL MIST ELIMINATOR KF25	Χ	949-9395
OIL MIST CARTRIDGE (0.ty 2)	X	949-9394

Tab. 14

RECOMMENDED OIL	P/N	Q.TY
Agilent Rotary Vane Fluid DS45 Type (Europe)	949-9390M001	1 liter

^(*) Contains all the valves, O-rings and seals to refurbish the pump to vacuum integrity.

^(**) Includes all the items of the minor maintenance kit plus the vanes.

Vacuum Products Division

Dear Customer,

Thank you for purchasing an Agilent vacuum product. At Agilent Vacuum Products Division we make every effort to ensure that you will be satisfied with the product and/or service you have purchased.

As part of our Continuous Improvement effort, we ask that you report to us any problem you may have had with the purchase or operation of our products. On the back side you find a Corrective Action request form that you may fill out in the first part and return to us.

This form is intended to supplement normal lines of communications and to resolve problems that existing systems are not addressing in an adequate or timely manner.

Upon receipt of your Corrective Action Request we will determine the Root Cause of the problem and take the necessary actions to eliminate it. You will be contacted by one of our employees who will review the problem with you and update you, with the second part of the same form, on our actions.

Your business is very important to us. Please, take the time and let us know how we can improve.

Sincerely

Giampaolo LEVI

Vice President and General Manager Agilent Vacuum Products Division

CUSTOMER REQUEST FOR CORRECTIVE / PREVENTIVE / IMPROVEMENT ACTION

AGILENT VACUUM PRODUCTS DIVISION TORINO – QUALITY ASSURANCE FAX

AGILENT TECHNOLOGIES ITALIA S.p.A. - Vacuum Products Division -

XXXX-011-9979350

TO:

ADDRESS:

Via F.Ili Varian, 54 – 10040 Leinì (TO) – Italy					
E-MAIL: vpd-qualityassurance_I	odl-ext@agilent.com				
NAME	COMPANY	FUNCTION			
ADDRESS:					
TEL. N°:	FAX N°:				
E-MAIL:					
PROBLEM / SUGGESTION :					
REFERENCE INFORMATION (model n	°, serial n°, ordering information, time	to failure after installation, etc.):			
NETERENCE IN ORWATION (Model II	, serial if , ordering information, time	to failure after installation, etc.,			
		ATE			
CORRECTIVE ACTION PLAN / ACTUAT (by AGILENT VPD)	ION	LOG N°			
(4)					

XXX = Code for dialing Italy from your country (es. 01139 from USA; 00139 from Japan, etc.)

Vacuum Products Division Instructions for returning products

Dear Customer,

Please follow these instructions whenever one of our products needs to be returned.

Complete the attached **Request for Return form** and send it to Agilent Technologies (see below), taking particular care to include the completed **Health and Safety** declaration Section. No work can be started on your unit until we receive a completed copy of this form.

After evaluating the information, Agilent Technologies will provide you with a **Return Authorization (RA) number** via email or fax, as requested. Note: Depending on the type of return, a Purchase Order may be required at the time **the Request for Return is submitted**. We will quote any necessary services (evaluation, repair, special cleaning, eg).

Product preparation

- Remove all accessories from the core product (e.g. inlet screens, vent valves).
- Prior to shipment and if applicable for your product, drain any oils or other liquids, purge or flush all gasses, and wipe off any excess residue.
- If ordering an Advance Exchange product, please use the packaging from the Advance Exchange to return the defective product.
- Seal the product in a plastic bag, and package product carefully to avoid damage in transit. You are responsible
 for loss or damage in transit.
- Include a copy of the Health and Safety Declaration in the shipping documentation on the outside of the shipping box of your returning product.
- Clearly label package with RA number. Using the shipping label provided will ensure the proper address and RA number are on the package. Packages shipped to Agilent without a RA clearly written on the outside cannot be accepted and will be returned.
- Return only products for which the RA was issued.

Shipping

- Ship to the location specified on the printable label, which will be sent, along with the RA number, as soon as
 we have received all of the required information. Customer is responsible for freight charges on returning
 product.
- Return shipments must comply with all applicable Shipping Regulations (IATA, DOT, ADR, etc.) and carrier requirements.

RETURN THE COMPLETED **REQUEST FOR RETURN** FORM TO YOUR NEAREST LOCATION: **EUROPE:**NORTH AMERICA: PACIFIC RIM:

Fax: 00 39 011 9979 330

Fax Free: 00 800 345 345 00 Fax: 1 781 860 9252 Please visit our website for Toll Free: 00 800 234 234 00 Toll Free: 800 882 7426, Option 3 individual office information

vpt-customercare@agilent.com vpl-ra@agilent.com http://www.agilent.com

Vacuum Products Division Terms and conditions

TERMS AND CONDITIONS

Please read the terms and conditions below as they apply to all returns and are in addition to the Agilent Technologies Vacuum Product Division – Products and Services Terms of Sale.

- Unless otherwise pre-negotiated, customer is responsible for the freight charges for the returning product.
 Return shipments must comply with all applicable Shipping Regulations (IATA, DOT, etc.) and carrier requirements.
- Agilent Technologies is not responsible for returning customer provided packaging or containers.
- Customers receiving an Advance Exchange product agree to return the defective, rebuildable part to Agilent Technologies within 15 business days. Failure to do so, or returning a non-rebuildable part (crashed), will result in an invoice for the non-returned/non-rebuildable part.
- Returns for credit toward the purchase of new or refurbished Products are subject to prior Agilent approval and may incur a restocking fee. Please reference the original purchase order number.
- O Units returned for evaluation will be evaluated, and a quote for repair will be issued. If you choose to have the unit repaired, the cost of the evaluation will be deducted from the final repair pricing. A Purchase Order for the final repair price should be issued within 3 weeks of quotation date. Units without a Purchase Order for repair will be returned to the customer, and the evaluation fee will be invoiced.
- Products returned that have not been drained from oil will be disposed.
- A Special Cleaning fee will apply to all exposed products.
- If requesting a calibration service, units must be functionally capable of being calibrated.

Vacuum Products Division Request for Return Form

Customer information					
Company:		Contact Name:			
Address:		Tel: Fax:		-	
Addices.		Email:		i ux.	
		Ciliali.			
Equipment			_		
Product description	Agilent PartNo	Agilent Serial No		Original Purchasing	Reference
Failure description		Type of process (fo	or which the equipr	nent was used)	
		,, , , , , , , , , , , , , , , , , , , ,			
Type of return					
☐ Non Billable ☐ Billable ➡ Ne	ew PO # (hard copy must be s	submitted with this form):	:		
,		Calibration		Return for Credit	
Health and safety		Substances (please	e refer to MSDS f	forms)	
·		· · ·		t that is exposed to radioact	ive, biological,
The product has been exposed to the fo	~			t written evidence of decont	-
(by selecting 'YES' you MUST complete the table to		Trade name	Chemical name	Chemical Symbol	CAS Number
Toxic	YES NO				
Harmful	YES NO		<u> </u>		
Corrosive	YES NO		<u> </u>		
Reactive	☐ YES ☐ NO				
Flammable	YES NO		<u> </u>		
Explosive (*)	YES NO		<u> </u>		
Radioactive (*)	YES NO		ļ		
Biological (*)	YES NO				
Oxidizing	YES NO				
Sensitizer	☐ YES ☐ NO				<u> </u>
Other dangerous substances	☐ YES ☐ NO				
Goods preparation					
If you have replied YES to one of the above questions of the above questions of the same o	uestions. Has the product b	een purged?		☐ YES ☐ NO	
Has the product been drained from oil?				YES NO	T APPLICABLE
I confirm to place this declaration on the outs	side of the shipping box.				I ALLEIONDEE
I declare that the ab	ove information is true	and complete to the	e best of my kno	wledge and belief.	
Lunderstand	d and agree to the term	es and conditions on	nage 2 of this d	locument	
T diffactorial in		- Contractions on	page 2 of tillo a	oddinent.	
Name:		Authorized Signature:			
Position:					
Date:					
NOTE: If a product is received at Agilent which is contaminated with a toxic or hazardous material that was not disclosed, the customer will be held responsible for all costs incurred to ensure the safe handling of the product, and is liable for any harm or injury to Agilent employees as well as to any third party occurring as a result of exposure to toxic or hazardous materials present in the product.					

Agilent Vacuum Products Division/Sales and Service Offices

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Agilent Technologies India Pvt. Ltd.

Unit Nos 105-116

First Floor, Splendor Forum, Plot No.-3 , District Centre, Jasola

New Delhi-110025 Ph: +91 11 4623 7100 Fax: +91 4623 7105 Toll Free: 18001801517

<u>Italy</u> Agilent Technologies Italia S.p.A.

Via F.lli Varian, 54
10040 Leini, (Torino) - Italy
Tel: +39 011 9979 111
Fax: +39 011 9979 350
Toll free: 00 800 234 234 00
vpt-customerservice@agilent.com

Janan

Agilent Technologies Japan, Ltd.

8th Floor Sumitomo Shibaura Building 4-16-36 Shibaura Minato-ku Tokyo 108-0023 - Japan Tel.: +81 3 5232 1253 Fax: +81 3 5232 1710

Toll-Free: 0120 655 040

vpj-customerservice@agilent.com

Korea

Agilent Technologies Korea, Ltd.

Shinsa 2nd Bldg. 2F, 966-5 Daechi-dong Kangnam-gu, Seoul Korea 135-280

Tel: +82 (0)2 2194 9449 Fax: +82 (0)2 3452 3947 Toll free: 080 222 2452

vpk-customerservice@agilent.com

Mexico

Agilent Technologies

Concepcion Beistegui No 109 Col Del Valle C.P. 03100 – Mexico, D.F.

Tel.: +52 5 523 9465 Fax: +52 5 523 9472

Singapore

Agilent Technologies Singapore Pte. Ltd,

1 Yishun Avenue 7, Singapore 768923 Tel: (65) 6215 8045 Fax: (65) 6754 0574 Toll free: 1 800 2762622

vps-customerservice@agilent.com

© Agilent Technologies, Inc. 2013 Printed in ITALY

12/2013

Publication Number: 87-900-143-01

Agilent Technologies

Southeast Asia

Agilent Technologies Sales Sdn Bhd

Unit 201, Level 2 uptown 2, 2 Jalan SS21/37, Damansara Uptown

47400 Petaling Jaya, Selangor, Malaysia Tel: +603 7712 6106 Fax: +603 6733 8121 Toll free: 1 800 880 805

vps-customerservice@agilent.com

Taiwan

Agilent Technologies Taiwan Limited

20 Kao-Shuang Road Ping-Chen City Tao-Yuan Hsien, 32450 Taiwan, R.O.C. Tel: +886 3 4959004 Toll free: 0800 018 768 vpw-customerservice@aqilent.com

UK and Ireland

Agilent Technologies UK, Ltd.

6 Mead Road Oxford Industrial Park Yarnton, Oxford OX5 1QU – UK Tel.: +44 (0) 1865 291570 Fax: +44 (0) 1865 291571

Other Countries

Agilent Technologies Italia S.p.A.

Via F.Ili Varian, 54 10040 Leini, (Torino) - Italy Tel.: +39 011 997 9111 Fax: +39 011 997 9350 Toll-Free: 00 800 234 234 00

Customer Support & Service

NORTH AMERICA:

Toll Free: 800 882 7426 vpl-ra@agilent.com

EUROPE:

Toll Free: 00 800 234 234 00 vpt-customercare@agilent.com

PACIFIC RIM:

please visit our website for individual office information

http://www.agilent.com/chem/vacuum

Worldwide Web Site, Catalog and Order

www.agilent.com/chem/vacuum Representatives in most countries

12/13